

## Gypsona

### Plaster of Paris Bandage

Gypsona bandages combine fast setting, high quality plaster of Paris with the conformable and stable interlocked weave of leno cloth. The bandage demonstrates minimum plaster loss whilst the rich creamy texture allows easy handling.

5cm	SN87060527	Box/12	\$48.48	2+	\$45.02
7.5cm	SN87067535	Box/24	\$104.13	2+	\$96.69
10cm	SN87061035	Box/24	\$134.79	2+	\$125.16
15cm	SN87061535	Box/12	\$94.49	2+	\$87.74



## Delta-Lite Plus

### Fibreglass Casting Bandage

Delta-Lite Plus is a fibreglass cast tape that combines a fibreglass substrate and advanced water activated polyurethane resin system which provides tack free properties. The tack free properties allow ease of application and molding using standard examination gloves.

- Easy to apply
- Greater conformability
- High lamination: smooth finished cast which is strong, long lasting and lightweight

5cm	SN7345801	Roll	\$10.07	10+	\$9.35
7.5cm	SN7345802	Roll	\$11.37	10+	\$10.56
10cm	SN7345803	Roll	\$13.85	10+	\$12.86
12.5cm	SN7345804	Roll	\$16.39	10+	\$15.22



Also available in following colours: red, blue, purple, black, pink, green, mixed  
Note that we only sell them in full cartons of 10. Price is equivalent to 10+ price of standard white.  
Colours are not a standard stock line

## Delta-Cast Elite

### Polyester Casting Tape

Delta-Cast Elite has been designed with ease of application and patient comfort in mind. Our patented fibreglass-free casting material combines strength and rigidity with soft edges for enhanced patient comfort.

- Smooth cast with soft edges
- Fibreglass free
- Easy to apply and remove
- Highly conformable

5cm	SN7227900	Roll	\$11.55	10+	\$10.73
7.5cm	SN7227901	Roll	\$13.08	10+	\$12.14
10cm	SN7227902	Roll	\$15.95	10+	\$14.81
12.5cm	SN7227903	Roll	\$19.38	10+	\$17.99



Also available in following colours:  
red, blue, purple, black, pink, green, mixed  
Note that we only sell them in full cartons of 10  
Price is equivalent to 10+ price of standard white.  
Colours are not a standard stock line



## Delta-Cast Print

### Polyester Casting Tape

Same features as Delta-Cast Elite but with an appealing cast finish

5cm	Camouflage	SN7227322	Box/10	\$178.78
5cm	Paw Print	SN7227320	Box/10	\$178.78
7.5cm	Camouflage	SN7227306	Box/10	\$198.06
7.5cm	Paw Print	SN7227321	Box/10	\$198.06



## 3M™ Scotchcast™ Plus

### Synthetic Casting Tape

Consists of a knitted fibreglass fabric impregnated with polyurethane resin. Exposure to moisture or water initiates a chemical reaction which causes the tape to become rigid. Finished casts produced are lightweight, strong, radiolucent, and porous. Scotchcast Plus casting tape contains a lubricating agent in the resin to provide a slippery feel when exposed to water.



5cm	3M82002	Roll	\$9.89	10+	\$9.18
7.5cm	3M82003	Roll	\$11.66	10+	\$10.83
10cm	3M82004	Roll	\$14.22	10+	\$13.20
12.5cm	3M82005	Roll	\$16.59	10+	\$15.41

## Delta-Cast Soft

### Semi-Rigid Polyester Casting Tape

Delta-Cast Soft is fibreglass-free and was designed to provide variable support to fractures and soft tissue injuries. Its semi-rigidity allows for more flexibility and mobility, thus helping to reduce immobilisation related problems such as atrophy or joint stiffness. As a result, the overall healing process tends to be shortened.

- Earlier patient mobility as cast will allow a limited range of movement
- Lightweight, breathable material offers high wearing comfort
- Smooth surface and edges for a snag-free finish

			White	Red	Blue
5cm	Box/10	\$173.56	SN7270812	SN7270808	SN7270805
7.5cm	Box/10	\$200.23	SN7270822	SN7270818	SN7270815
10cm	Box/10	\$234.72	SN7270832	SN7270828	SN7270825



BSN medical

## Protouch

### Synthetic and Cotton Stockinette

- Soft and comfortable providing good cushioning and comfort to the patient.
- Stretches to three times its normal width to conform to an arm, leg or body.

#### Protouch Synthetic Stockinette

- Made of synthetic fibres to react moisture away from the skin for rapid evaporation

5cm	SN301002	Roll	\$28.03
7.5cm	SN301003	Roll	\$33.73
10cm	SN301004	Roll	\$57.01



#### Protouch Cotton Stockinette

- Ideal material for an initial undercast layer
- Made of a high-grade knitted cotton

5cm	SN307002	Roll	\$22.22
7.5cm	SN307003	Roll	\$25.48
10cm	SN307004	Roll	\$30.64



## Tubular-Rib™

### Tubular Retention Bandage

Designed to provide a protective layer in direct contact with the skin

#### Used for:

- Protection of the skin prior to the application of fibreglass or plaster of paris casts.
- Protection of the skin under all types of compression bandaging.

#### Features:

- Australian made
- 100% Cotton

5cm	Small Limbs	SUTSM500	\$17.64
7.5cm	Medium Limbs	SUTSM501	\$23.80
10cm	Large Limbs	SUTSM502	\$35.56
15cm	Extra Large Limbs	SUTSM503	\$50.54

Size is unstretched width  
All rolls are 20m



#### Application:

It is important when measuring Tubular-Rib prior to application that an amount of 5cm is allowed at either end of the bandage. The bandage can be gently rolled onto the limb or torso allowing a split at the thumb site for arm casts.

Tubular-Rib can be applied independently using multiple layers or in conjunction with flat undercast padding.

## Tubular-Plus™

### Tubular Padded Protection Bandage

Designed to provide a protective layer in direct contact with the skin.

#### Used for:

- Protection of the skin prior to the application of fibreglass & plaster of paris casts.
- Padding to limbs under compression bandaging.

#### Features:

- Australian made
- Cotton 98% (Elastodiene / Nylon 2%)

5cm	Small Limbs	SUTSM6629	\$42.65
7.5cm	Medium Limbs	SUTSM6630	\$45.24
10cm	Large Limbs	SUTSM6631	\$56.60
12.5cm	Extra Large Limbs	SUTSM6632	\$59.50

Size is unstretched width  
All rolls are 5m



#### Application:

It is important when measuring Tubular-Plus prior to application, that an amount of 2.5cm is allowed at either end of the bandage. The 2.5cm of additional bandage will gently roll after being cut to prevent airborne particles from entering the wound.

After using Tubular-Plus for an extended period, the bandage will be subtly shaped to the limb it has been applied to. Gently pull the bandage lengthways to remove the memory and the bandage will return to its original specifications.

## Soffban Natural

### Orthopaedic Padding

- Conforms well for easy shaping
- Natural fleece absorbs sweat and exudate
- Ideal for plaster casts
- 2.7m



5cm	SN7147200	Pkt/12	\$24.35	2+	\$22.61
7.5cm	SN7147300	Pkt/12	\$25.83	2+	\$23.99
10cm	SN7146900	Pkt/12	\$32.91	2+	\$30.56
15cm	SN7147000	Pkt/12	\$38.44	2+	\$35.70

## Soffban Synthetic

### Water-Resistant Orthopaedic Padding

- Water-resistant
- Ideal for synthetic casts
- 2.7m



5cm	SN36000025	Pkt/12	\$24.35	2+	\$22.61
7.5cm	SN36000026	Pkt/12	\$25.83	2+	\$23.99
10cm	SN36000027	Pkt/12	\$32.91	2+	\$30.56
15cm	SN36000028	Pkt/12	\$38.44	2+	\$35.70

## Velband

### Absorbent Orthopaedic Padding (2.7m)

Velband is a hypoallergenic underbandage with a unique three-layer structure for easy application and long-term comfort.

It provides excellent cushioning and moisture management under plaster of Paris or synthetic casts, avoiding irritation and pressure sores by keeping the skin dry and comfortable.



7.5cm	SN7231501	Pkt/12	\$24.25	2+	\$22.52
10cm	SN7231502	Pkt/12	\$29.16	2+	\$27.08
15cm	SN7231504	Pkt/12	\$42.98	2+	\$39.91



## Delta-Dry

### Waterproof Cast Liner

- Water resistant cast padding (2.4m)
- Allows patients to shower, swim and bathe like normal
- Better long term water resistance
- Offers more under cast protection
- Eliminates the need for stockinette



5cm	SN7344300	Box/12	\$219.20
7.5cm	SN7344301	Box/12	\$264.14
10cm	SN7344302	Box/12	\$309.19

## Pure Cotton

### Orthopaedic under-bandages

Made from high quality, long staple cotton fibres.

Designed and made to have some self-adhesion property and appropriate tensile strength, resistant to abrasion and plaster strike through.

2.7m

5cm	MGT30700	Pkt/12	\$10.95	3+	\$9.82
7.5cm	MGT30701	Pkt/12	\$15.42	3+	\$13.84
10cm	MGT30702	Pkt/12	\$16.24	3+	\$14.58
15cm	MGT30703	Pkt/12	\$26.81	3+	\$24.06



## Cotton undercast padding

Ideal under plaster casts - 2.7m

5cm	AXP10105731	Pkt/12	\$8.34	3+	\$6.99
7.5cm	AXP10105732	Pkt/12	\$9.37	3+	\$8.44
10cm	AXP10105733	Pkt/12	\$11.46	3+	\$10.50
15cm	AXP10105734	Pkt/6	\$12.60	3+	\$11.70

BUDGET SAVER



## Viscose undercast padding

Ideal under synthetic casts - 2.7m

5cm	AXP10105711	Pkt/12	\$8.02	3+	\$6.00
7.5cm	AXP10105712	Pkt/12	\$9.04	3+	\$8.06
10cm	AXP10105713	Pkt/12	\$10.08	3+	\$9.12
15cm	AXP10105714	Pkt/6	\$11.02	3+	\$8.25



## Finger Stax splints

Sz1	MIAS01	Sz5	MIAS05
Sz2	MIAS02	Sz6	MIAS06
Sz3	MIAS03	Sz7	MIAS07
Sz4	MIAS04		
Each	\$4.25		

## Finger splints

<b>Aluminium</b>			
13mm x 480mm	ROW4813	\$3.45	2+ \$3.22
19mm x 480mm	ROW4819	\$3.75	2+ \$3.50
24mm x 480mm	ROW4824	\$4.05	2+ \$3.78
32mm x 480mm	ROW4832	\$4.50	2+ \$4.20

## Sam Splint

Lightweight, reusable and very portable

The Sam Splint Original is an extremely pliable multi-purpose splint and is built from a thin core of aluminium alloy, sandwiched between two layers of closed cell foam. Bent into various shapes and curves, it becomes extremely strong and supportive for any fractured or injured limb.

9.5cm x 4cm	FERSVPSSF	\$2.95
91cm x 10mm	FERSVPSS36	\$19.65



**FERNO**  
Australia

## Ferno Vacuum Splint Set

Each splint is constructed of durable vinyl coated nylon that is supple enough to be moulded around an injured extremity.

Air is extracted via the pump supplied forcing the splint to become hard and mould around the patient providing local immobilisation without inward pressure, reducing the risk of additional trauma.

The complete set includes three splints (small, medium and large), compact vacuum pump, repair kit and carry case.

FERSVVSF  
**\$949**

Provides maximum support with minimum effort



**FERNO**  
Australia

## Dynacast Prelude Precuts

Acute Synthetic Splinting System

Dynacast Prelude is an easy-to-use, pre-assembled synthetic splinting system. Made from polyurethane coated fibreglass covered by polypropylene padding, Dynacast Prelude provides stabilisation in 2-4 minutes and enables weight bearing in just 20 minutes

5 x 25cm	SN7286900	Box/5	\$52.49
7.5 x 30cm	SN4722200	Box/5	\$85.41
10cm x 40cm	SN4722500	Box/5	\$135.92
12.5 x 76cm	SN4723500	Box/4	\$257.59



**BSN** medical

## Delta-Xpress

Prefabricated Splinting System

Delta-Xpress is a groundbreaking new splinting system that allows you to immobilise limbs within a fraction of the time when compared to conventional splints or casts. The prefabricated design is intuitive to use and eliminates many time-consuming steps while still providing secure and strong immobilisation.

- Time saver - fast immobilisation
- Easy to apply
- Immobilisation - secure, safe and strong

### Wrist Splint

Medium	SN7271201	Each	\$16.43	2+	\$15.33
Large	SN7271202	Each	\$16.43	2+	\$15.33

### Elbow

Small	SN7271000	Each	\$19.74	2+	\$18.42
Medium	SN7271001	Each	\$19.74	2+	\$18.42
Large	SN7271002	Each	\$19.74	2+	\$18.42



**BSN** medical



## Slings

Disposable Slings 110 x 110 x 155cm	AXP10106311	Each	\$0.49	12+	\$0.40
Cotton Slings 110 x 110 x 155cm	AXP10106211	Each	\$0.87	12+	\$0.69
Safety Pins 57mm	AXP10110414	Pkt/12	\$0.51		

## Collar & Cuff Arm Sling



6m  
ROWCC06W  
\$12.15

## Actimove Sling

### Collar and Cuff Support

The padded strap and Y-tab closure system provide effective vertical support to the collarbone and upper limb area

- Extra patient comfort with padded strap
- Y-tab hook system facilitates secure clinical positioning
- Convenient dispenser box presentation
- 28 Y-tabs are included in each box

5cm x 12m SN7285918 Box/2rolls \$43.60



## Actimove Mitella

### Comfort Arm Sling

The envelope sling and neck strap of the Actimove Mitella provide effective vertical support to the arm.

- Extra patient comfort with neck pad
- Adjustable neck strap, with Y-tab hook system facilitates secure clinical positioning
- Soft thumb loop maintains sling placement

Small SN7281961  
Medium SN7281962  
Large SN7281963  
X-Large SN7281964 \$50.06



## Actimove Umerus

### Comfort Shoulder Immobiliser

Actimove Umerus is a shoulder immobiliser that limits abduction, extension and internal rotation of the humerus

- Integrated neck pad provides additional patient comfort
- Adjustable neck and waist strap with Y-tab hook system, facilitates secure clinical positioning
- Elbow mobilisation is possible without complete sling removal
- Soft thumb loop effectively maintains sling placement

Small SN7281939 Large SN7281941  
Medium SN7281940 X-Large SN7281942 \$58.62



## Actimove Genu Tri Panel

### Three Panel Knee Immobiliser

Provides effective immobilisation of the knee and helps to reduce pain following trauma or surgery. The 'tri-panel' design allows splint adjustments for a broad spectrum of patient leg circumferences and the hook and loop sensitive closure system permits easy application and removal.

50cm Small SN7281914 \$77.22  
50cm Large SN7281915 \$83.05  
60cm Small SN7281916 \$88.91  
60cm Large SN7281917 \$95.58



## Actimove Manus Forte

### Wrist Immobiliser

A highly breathable, stabilising wrist brace for the treatment of painful conditions and irritations of the wrist joint.

Small/Med Right SN7348203  
Small/Med Left SN7348204  
Lrg/X-Lrg Right SN7348205  
Lrg/X-Lrg Left SN7348206

\$33.98



## Cast Shoe

Constructed from a durable upper canvas and cushioned foam. Combined with pressure sensitive closures and slip-resistant, rocker shaped sole, these shoes offer comfort and security to the wearer.

- Helps promote normal gait pattern
- Three hook and loop closures ensure a precise fit

X-Small 19 x 10.5cm SN7981112 \$31.36  
Small 22 x 10.5cm SN7981113 \$31.36  
Medium 26.5 x 11.5cm SN7981115 \$37.67  
Large 30 x 13cm SN7981117 \$37.67



## Cervical collar

### Ferno Wizloc collar

Universal design provides a collar capable of fitting 90% of patients. With its three step adjustment, application is simple and fast.

- Designed for use on 2 year olds through to adult
- On-patient sizing for greater accuracy
- Can be decontaminated and reused up to 10 times
- Three position occipital support with keyhole lock pin
- X-ray transparent and MRI compatible



FERCCFWLG

\$42\*

FERCCFWL

## Cervical Collar

Small AXP10029501  
Medium AXP10029502  
Large AXP10029503

\$7.63\*

