

## SurGuard2® Safety Needles

Incorporates a simple locking shield which significantly reduces the risk of needle stick injuries. The activation of the safety device can be easily performed with one hand.

Box/100 \$21.55\*

**safety**

- 23G x 1" TER6SG22325
- 23G x 1.5" TER6SG22338
- 25G x 1" TER6SG22525




## VanishPoint® Retractable Syringes

- Reduces exposure - Passive retraction, needle retracts directly from patient into syringe
- Easy to use single-handed

23G x 1" TER8RTI10311  
 25G x 1" TER8RTI10391  
 30G x 1/2" TER8RTI15221

Box/100 \$55.13\*

**safety**



## Surshield Safety Winged Infusion Sets

- Revolutionary design in sharps injury protection
- Protector is snapped over the needle immediately after its withdrawal from patient's vein
- Single handed technique

9cm	19G	TER5SVS19BLS	30cm	19G	TER5SVS19BL
9cm	21G	TER5SVS21BLS	30cm	21G	TER5SVS21BL
9cm	23G	TER5SVS23BLS	30cm	23G	TER5SVS23BL
9cm	25G	TER5SVS25BLS	30cm	25G	TER5SVS25BL

**safety**

Box/50 \$70.29\*



## Surflo Winged Infusion Sets Scalpvein Sets

9cm	19G	TER1SV19BLS	30cm	19G	TER1SV19BLK
9cm	21G	TER1SV21BLS	30cm	21G	TER1SV21BLK
9cm	23G	TER1SV23BLS	30cm	23G	TER1SV23BLK
9cm	25G	TER1SV25BLS	30cm	25G	TER1SV25BLK

Each \$0.84\*      100+ \$0.75\*



## Surflo IV Catheters

The tip of the introducing needle is cut and honed to an ultra-sharp triple-bevel that is designed to minimise patient discomfort and tissue trauma. The gently tapered catheter aligns with the needle bevel to provide smooth transition from needle to catheter during insertion.

14G x 2"	TER6SROX1451C	20G x 1.25"	TER6SROX2032C
16G x 2"	TER6SROX1651C	20G x 2"	TER6SROX2051C
18G x 1.25"	TER6SROX1832C	22G x 1"	TER6SROX2225C
18G x 2"	TER6SROX1851C	24G x 3/4"	TER6SROX2419C
20G x 1"	TER6SROX2025C		

Each \$1.34\*      50+ \$1.24\*      100+ \$1.19\*





## Syringes

1ml	Slip Tuberculin Tip	TER6SS01T	Box/100	\$16.23*	5+	\$15.05*
3ml	Slip Tip	TER6SS03S	Box/100	\$9.73*	5+	\$9.02*
5ml	Slip Tip	TER6SS05S	Box/100	\$14.45*	5+	\$13.40*
10ml	Slip Tip	TER6SS10S	Box/100	\$18.69*	5+	\$17.32*
3ml	Lock Tip	TER6SS03L	Box/100	\$11.59*	5+	\$10.74*
5ml	Lock Tip	TER6SS05L	Box/100	\$17.15*	5+	\$15.90*
10ml	Lock Tip	TER6SS10L	Box/100	\$22.88*	5+	\$21.21*
20ml	Lock Tip	TER6SS20L	Box/50	\$24.84*	5+	\$23.03*
30ml	Lock Tip	TER1SS30LE1	Box/50	\$41.58*	5+	\$38.54*
50ml	Lock Tip	TER1SS50LE	Box/20	\$24.65*	5+	\$22.85*
20ml	Slip Eccentric Tip	TER6SS20ES	Box/50	\$17.89*	5+	\$16.59*
30ml	Slip Eccentric Tip	TER1SS30ESE1	Box/50	\$43.40*	5+	\$40.23*
50ml	Slip Eccentric Tip	TER1SS50ESE	Box/20	\$25.55*	5+	\$23.69*
50ml	Catheter Tip	TER1SS50CE	Box/20	\$27.01*	5+	\$24.34*



## Needles

18G x 1.5"	TER1NN1838R	23G x 1.25"	TER1NN2332R
19G x 1.5"	TER1NN1938R	25G x 1"	TER1NN2525RL
21G x 1"	TER1NN2125RL	25G x 3/4"	TER1NN2519RA
21G x 1.25"	TER1NN2132R	25G x 5/8"	TER1NN2516RA
21G x 1.5"	TER1NN2138R	26G x 1/2"	TER1NN2613R
22G x 1.5"	TER1NN2238R	27G x 1/2"	TER1NN2713R
23G x 1"	TER1NN2325R	30G x 1/2"	TER3NN3013

Box/100 \$4.80\* 10+ \$4.45\*



### Drawing Up Needle 18G x 1.5"

TER1XXNN1838  
Box/100 \$5.08\*

## Syringes with needles

3ml 23G x 1.25"	TER6SS03S2332	Box/100	\$21.59*	5+	\$20.02*
3ml 25G x 5/8"	TER6SS03S2516	Box/100	\$21.59*	5+	\$20.02*
5ml 21G x 1.5"	TER6SS03S2138	Box/100	\$32.39*	5+	\$30.02*
5ml 23G x 1.25"	TER6SS05S2332	Box/100	\$26.39*	5+	\$24.46*



## Insulin syringes

0.5ml 27G x 1/2"	TER3SS05M2713A	Box/100	\$21.59*	5+	\$20.02*
0.5ml 29G x 1/2"	TER3SS05M2913A	Box/100	\$22.45*	5+	\$20.82*
1ml 27G x 1/2"	TER3SS10M2713A	Box/100	\$21.59*	5+	\$20.02*
1ml 29G x 1/2"	TER3SS10M2913A	Box/100	\$22.45*	5+	\$20.82*
1ml Slip Tip	TER6SS01NA	Box/100	\$15.44*	5+	\$14.31*



Terumo insulin syringes are designed to provide minimal patient discomfort, accurate dosing and consistent high quality

## Spinal needles



18G x 90mm	TER1SN1890
19G x 90mm	TER1SN1990
20G x 38mm	TER1SN2038
20G x 70mm	TER1SN2070
20G x 90mm	TER1SN2090
21G x 70mm	TER1SN2170
21G x 90mm	TER1SN2190
22G x 38mm	TER1SN2238
22G x 70mm	TER1SN2270
22G x 90mm	TER1SN2290
23G x 70mm	TER1SN2370
23G x 90mm	TER1SN2390
25G x 90mm	TER1SN2590

38mm = 1 1/2"  
70mm = 2 3/4"  
90mm = 3 1/2"

Box/20 \$40.83\*  
2+ \$37.85\*

## Dental needles

Terumo dental needles have an ultra sharp lancet point designed to minimise pain and discomfort

25G x 21mm	TER1DN2521
25G x 41mm	TER1DN2541
27G x 22mm	TER1DN2722
27G x 35mm	TER1DN2735
27G x 41mm	TER1DN2741
30G x 22mm	TER1DN3022

Box/100 \$12.23\*  
5+ \$11.34\*



## Eclipse needles

Safety within your reach

- Shielding safety technology
- Activates with a single finger stroke immediately after injection, no hard surface required
- For use with BD Luer-Lock syringes only
- ^ SmartSlip™ Technology (ok with luer slip)

<b>18g x 1.5"</b>	BD305766	<b>25g x 1"</b>	BD305761
<b>21g x 1"</b>	BD305764	<b>25g x 1.5"</b>	BD305767
<b>21g x 1.5"</b>	BD305765		
<b>22g x 1"</b>	BD305768	<b>25g x 5/8"</b>	BD305760 ^
<b>22g x 1.5"</b>	BD305763	<b>27g x 1/2"</b>	BD305570 ^
<b>23g x 1"</b>	BD305762	<b>30g x 1/2"</b>	BD305771 ^



Box/100 \$33.57\*  
5+ \$31.12\*

## Insyte IV Catheters



<b>14G x 45mm</b>	BD381267		
<b>16G x 30mm</b>	BD381254		
<b>16G x 45mm</b>	BD381257		
<b>18G x 30mm</b>	BD381244		
<b>18G x 48mm</b>	BD381247		
<b>20G x 30mm</b>	BD381234		
<b>20G x 48mm</b>	BD381237		
<b>22G x 25mm</b>	BD381223	Each	<b>\$1.98*</b>
<b>24G x 19mm</b>	BD381212	50+	<b>\$1.67*</b>

## Insyte Autoguard Shielded IV Catheters

Designed to reduce needlestick injuries



The unique push-button shielding mechanism releases the spring and allows the needle and flash chamber to quickly retract into the safety barrel. The clinician maintains control of the process by deciding when to activate the push-button shielding mechanism.

<b>16G x 30mm</b>	BD381454		
<b>18G x 30mm</b>	BD381444		
<b>18G x 48mm</b>	BD381447		
<b>20G x 25mm</b>	BD381433		
<b>20G x 30mm</b>	BD381434	Box/50	<b>\$142.50*</b>
<b>22G x 25mm</b>	BD381423	2+	<b>\$136.80*</b>

## BD Insyte Autoguard BC

With Blood Control technology

Blood exposure during a peripheral IV insertion is still a risk factor to healthcare workers and patients

- Reduces blood leakage: less need to worry about stopping blood
- Less cross-contamination potential
- Less risk of exposure to bloodborne pathogens
- Reduced need to clean up blood spills
- Cut your clean-up supply costs

## Making safety safer



BD Blood Control Technology - stops blood so it doesn't escape the hub. Controlling blood escape means a cleaner area helping to increase patient confidence.

<b>16G x 30mm</b>	BD382554	<b>18G x 48mm</b>	BD382547	<b>20G x 48mm</b>	BD382537		
<b>16G x 45mm</b>	BD382557	<b>20G x 25mm</b>	BD382533	<b>22G x 25mm</b>	BD382523	Box/50	<b>\$152.07*</b>
<b>18G x 30mm</b>	BD382544	<b>20G x 30mm</b>	BD382534	<b>24G x 19mm</b>	BD382512	2+	<b>\$145.99*</b>

## BD PosiFlush SP

Pre-filled Saline Syringe

These pre-filled syringes contain sodium chloride 0.9% and are designed for flushing in-situ vascular access devices. The 3ml, 5ml and 10ml syringes all have the same diameter as a standard 10ml syringe. May be used for flushing catheters when protocols stipulate use of a 10ml syringe to minimise flush pressure (PSI).

- Luer-Lok™ tip for secure connection to the catheter
- Eliminates syringe-induced blood reflux, enhancing catheter maintenance protocol
- Size and shape of the tip cap allow easy handling, ensures no touch contamination of the syringe luer connection
- Eliminates the use of needles, thus reducing the potential for accidental needlestick injuries
- Standardised flush volumes
- Must not be used on a sterile field

<b>3ml</b>	BD306573	Box/30	<b>\$14.36</b>
<b>5ml</b>	BD306574	Box/30	<b>\$15.57</b>
<b>10ml</b>	BD306575	Box/30	<b>\$17.28</b>





## Syringes

1ml	Slip Tuberculin Tip	BD302100	Box/100	\$14.11*	8+	\$13.08*
2ml	Slip Tip	BD302204	Box/100	\$9.32*	8+	\$8.64*
5ml	Slip Tip	BD302130	Box/100	\$14.11*	4+	\$13.08*
10ml	Slip Tip	BD302143	Box/100	\$17.37*	4+	\$16.10*
20ml	Slip Tip	BD300140	Box/50	\$22.67*	4+	\$21.02*
1ml	Lock Tip	BD309628	Box/100	\$54.66*	5+	\$50.67*
3ml	Lock Tip	BD302113	Box/100	\$10.97*	4+	\$10.17*
5ml	Lock Tip	BD302135	Box/100	\$15.89*	4+	\$14.73*
10ml	Lock Tip	BD302149	Box/100	\$20.56*	4+	\$19.06*
50ml	Lock Tip	BD300144	Box/40	\$45.43*	2+	\$42.11*
10ml	Slip Eccentric Tip	BD302146	Box/100	\$17.40*	4+	\$16.13*
20ml	Slip Eccentric Tip	BD300613	Box/120	\$53.10*	4+	\$49.23*
30ml	Slip Eccentric Tip	BD301231	Box/60	\$37.81*	4+	\$35.05*
50ml	Slip Eccentric Tip	BD300145	Box/40	\$41.01*	4+	\$38.02*
50ml	Catheter Tip	BD300867	Box/60	\$65.27*	4+	\$60.50*



## Syringes with needles

3ml	23G x 1.25"	BD302119	Box/100	\$21.59*	5+	\$20.02*
5ml	23G x 1.25"	BD302137	Box/100	\$30.69*	5+	\$28.45*



## Needles

18G x 1.5"	BD302032	23G x 1.25"	BD302008		
19G x 1.5"	BD301828	25G x 5/8"	BD301805		
20G x 1"	BD302025	25G x 1"	BD301807		
20G x 1.5"	BD302026	25G x 1.5"	BD301808		
21G x 1.5"	BD302017	26G x 1/2"	BD302002	Box/100	\$4.55*
22G x 3/4"	BD301815	26G x 5/8"	BD304300	10+	\$4.16*
23G x 3/4"	BD301810	27G x 1/2"	BD301801		
23G x 1"	BD302006	30G x 1/2"	BD304000		



## Other needles

Needle 16G x 1.5"	BD300637
Box/100	\$7.95*
10+	\$6.94*

Drawing Up Needles 18G x 1.5"	BD300204
Box/100	\$5.38*
10+	\$4.84*

### Blunt Fill Needle 18G x 1.5"

- Specifically designed for safe injection preparation
- Eliminates an unnecessary sharp
- 10 x the force to penetrate the skin, compared to a conventional needle

Box/100	\$23.36*	
10+	\$21.63*	BD305180



## Ultra-Fine Insulin syringes



0.3ml	29G x 1/2"	BD326103	Box/100	\$24.39*	5+	\$22.61*
0.3ml	31G x 5/16"	BD328822	Box/100	\$24.39*	5+	\$22.61*
0.5ml	29G x 1/2"	BD326769	Box/100	\$25.21*	5+	\$23.37*
0.5ml	30G x 1/2"	BD326725	Box/100	\$25.21*	5+	\$23.37*
0.5ml	31G x 5/16"	BD328821	Box/100	\$25.10*	5+	\$23.32*
1ml	29G x 1/2"	BD326719	Box/100	\$25.21*	5+	\$23.37*
1ml	31G x 5/16"	BD328820	Box/100	\$24.39*	5+	\$22.61*

## IV Catheters

See pages 256 & 258



## Giving Sets



### Straight Type Blood/Solution Set with 200 Micron Filter, Interlink

Blood/solution set, single non-vented drip chamber, 200 micron filter, roller clamp, 1 Interlink injection site 15cm from luer lock adapter (capped) with a retractable collar.



**BXTFNC2211**  
 Each \$6.55  
 48+ \$6.10

### Blood/solution set long line hand pump, Interlink

Blood Solution Set, non-vented drip chamber, 200 micron filter, hand pressure pump, roller clamp, 2 Interlink injection sites 30cm and 200cm from luer lock adapter (capped) with a retractable collar.



**BXTFNC1144**  
 Each \$23.00  
 50+ \$21.41

### Continue-Flo Solution Set, Clearlink

Solution Continue-Flo set, non-vented drip chamber with 15 micron filter, non-return valve, roller clamp, 2 Clearlink luer activated valves 15cm and 200cm from luer lock adapter (capped) with a retractable collar



**BXTFNC3120**  
 Each \$10.28  
 50+ \$9.57



## Giving sets



Transfusion set, 220cm length, vented 20 drops/ml chamber, roller regulator, rotating male luer lock collar, swabable "Robertsite" needle free injection site

**MMD020001**  
 Each \$5.45      50+ \$4.85

## Injection site

Interlink injection site with luer lock adapter

**BXT2N3399K**  
 Each \$2.36  
 10+ \$2.02



## Blunt Plastic Cannula

A Truly Universal Cannula

Designed for improved performance, the Blunt Plastic Cannula enhances the Interlink® IV Access System. The Blunt Plastic Cannula is a truly universal cannula, and can be used with all types of split septum IV injection sites and vials designed for needleless access.

**BD303345**  
 Box/100 \$34.25  
 2+ \$31.50



## Q-Syte™ Luer Access Split Septum

Ultimate performance in closed luer access  
 A simple split-septum penetrated by a standard male luer (syringe tip or administration set) to become the fluid path

**BD385100**  
 Box/50 \$145.01



## Stopcock 3 Way with Luer Lock

**BD394600**  
 Each \$1.16  
 50+ \$0.96

