

Breathe better with **BREATH-A-TECH®**



Upgraded **Autoclavable Hospital Grade Spacer**



**Low resistance airflow via
Patented Cross Valve
Technology.™**



Robust & fully autoclavable.



Clinically tested efficacy.



Designed for hospital applications.



Latex & BPA free.



Suitable for multi patient use.



Proudly Australian designed & owned.



SCAN FOR MORE INFO

Customer Information (Toll Free in Australia): 1800 241 596
Email: info@breathatech.com
Explore the full range at www.breathatech.com



**BREATH
A-TECH®**

The Breath-A-Tech® Hospital Grade Spacer is a hand-held, portable valved holding chamber (VHC). The spacer is placed between a pressurised metered dose inhaler (pMDI, or "puffer") and patient's mouth to function as a reservoir from which aerosol medication is dispensed. In order to minimise the delivery of large, aerosolised particles, the spacer reduces the direct delivery of large, aerosolised particles into the mouth, which are intended for the lower airways. Suitable for adults and children to use with most standard puffers in the treatment of asthma and COPD. ALWAYS FOLLOW THE DIRECTIONS FOR USE.