

gloveon Matador

Nitrile Exam Gloves Powder Free, Long Cuff



GloveOn® Matador is built with quality, consistency and reliability in an 8ml nitrile glove. With GloveOn® Matador, users can take confidence in wearing a stronger glove to shape their performance in their work duties.



GloveOn® Matador		
Length (mm)	≥ 280	
Thickness Measurements (mm)		
Palm (centre of Palm)	0.20 ± 0.04	
Finger (13mm ± 3mm from tip)	0.28 ± 0.04	
Cuff (25mm from cuff edge)	0.14 ± 0.04	
Physical Properties	Before Ageing	After Ageing
Tensile Strength (MPa)	≥ 18	≥ 16
Elongation (%)	≥ 500	≥ 400
Inspection Levels & AQL	Inspection Level	AQL
Watertightness	G1	1.50
Physical Dimensions	S2	4.00
Tensile Strength	S2	4.00
Visual Inspection (Major)	S4	2.50
Visual Inspection (Minor)	S4	4.00
Particulate Residue	N = 5	≤ 2mg/glove

Tested in accordance with EBS EN 16523-1:2016 with the following chemical permeation performance levels:

Chemotherapy Drugs and Concentration (Tested for Resistance to Permeation by Chemotherapy Drugs as per ASTM D6978-05 Test Report PN 107200B)	Minimum Breakthrough Detection Time (minutes)	EN ISO 374-1:2016 Type A	Level
Carmustine (BCNU), 3.3mg/ml (3,300 ppm)	90.6 minutes	Sulphuric Acid 96% (L)	2
Cisplatin, 1.0mg/ml (1,000 ppm)	>240 minutes	n-Heptane (J)	3
Cyclophosphamide (Cytoxan), 20.0mg/ml (20,000 ppm)	>240 minutes	Sodium Hydroxide 40% (K)	6
Dacarbazine (DTIC), 10.0mg/ml (10,000 ppm)	>240 minutes	Ammonium hydroxide 25% (O)	4
Doxorubicin Hydrochloride, 2.0mg/ml (2,000 ppm)	>240 minutes	Formaldehyde 37% (T)	6
Etoposide (Toposar), 20.0mg/ml (20,000 ppm)	>240 minutes	Hydrogen peroxide 30% (P)	6
Fluorouracil, 50.0mg/ml (50,000 ppm)	>240 minutes	Hydrofluoric acid 40% (S)	2
Methotrexate, 25.0mg/ml (25,000 ppm)	>240 minutes		
Mitomycin C, 0.5mg/ml (500 ppm)	>240 minutes		
Paclitaxel (Taxol), 6.0mg/ml (6,000 ppm)	>240 minutes		
Thiotepa, 10.0mg/ml (10,000 ppm)	165.8 minutes		
Vincristine Sulfate, 1.0mg/ml (1,000 ppm)	>240 minutes		

WARNING: Carmustine and Thiotepa, at the tested concentration, degraded Matador nitrile glove at 90.6 minutes and 165.8 minutes, respectively. The safe use of gloves in chemotherapy treatment is solely the decision of clinicians authorised to make such decision.

REORDER CODE

NTR18SS	SMALL
NTR18MM	MEDIUM
NTR18LL	LARGE
NTR18XL	X-LARGE
NTR18XXL	XX-LARGE

FEATURES

- Fingertip textured
- Powder free
- Not made with natural rubber latex
- Chemo drugs tested
- Lab chemical tested
- Ambidextrous
- Long cuff
- Sky blue colour

PACKAGING

50 gloves per box for S to XL
45 gloves per box for XXL
10 boxes per carton

REGULATORY COMPLIANCE

ARTG 164563, FDA 510(k),
MDD 93/42/EEC, REACH,
ROHS Directive 2002/95/EC,
EU 2016/425

STANDARDS

AS/NZS 4011-1, ASTM D6319, ASTM D5151,
ASTM D6124, ASTM D6978, EN 420,
EN 16523-1, EN ISO 374-1 (Type A),
EN 374 part 2, 4 & 5, ISO 11193-1,
EN 455 part 1, 2, 3 & 4, ISO 10993 part 5 & 10,
HACCP Certified

MANUFACTURING ACCREDITATIONS

ISO 9001, ISO 13485

A brand by



Mun Australia
Toll free: 1800 456 837
munglobal.com.au



Proud Supporter of
Birthing Kit Foundation (Australia)



Protection Always On