



Medisafe Biomedical 258ULT



Why Choose us

Medisafe Plus is the original premium brand of vaccine fridge within Australia. It is designed and built in Australia and has been protecting vaccines and medical products since 2010.

Medisafe Biomedical storage units are purpose-built for use in hospitals and laboratories. The units are designed to ensure low noise and consistent temperature readings through a fridge and freezer range.

A Medisafe Biomedical storage unit is feature packed, resulting in a storage solution that gives you peace of mind - knowing your vaccines/drugs/samples are stored in the safest possible manner.

Medisafe Biomedical can tailor our solutions to end user needs to allow for additional probes or temperature recording systems through the manufacturing process.

Medisafe Plus cloud monitoring is optional on all Medisafe Biomedical units, with 12 months free access to the portal. The cloud monitoring acts as an automated data logger which allows you to view and print your temperature chart with Min and Max and data in real time on any web enabled smartphone, tablet, or computer from any location.

All units have a BMS system accessible in the back of the unit for hospital remote monitoring.

Medisafe Biomedical 258ULT

Medisafe Biomedical 258ULT

Is an upright ultra low temperature unit in the Medisafe Biomedical range.

Medisafe Biomedical 258ULT medical freezer is a purpose-built medical freezer designed for the holding of samples and vaccines that require ultra low temp holding

Medisafe Biomedical fridges/freezers are feature packed, resulting in a storage solution that gives you peace of mind – knowing your vaccines/drugs/samples are stored in the safest possible manner.

The Medisafe Biomedical 258ULT is TGA approved.

The optional cloud monitoring acts as an automated data logger which allows you to view and print your temperature data in real time, from any location, on any web enabled smartphone, tablet, or computer.

All units have a BMS system plugged into the back of the unit for hospital remote monitoring.



Features

- USB Connection
- Temperature Graphs
- Temperature Controller
- Audio/Visual/Probe/ Open Door/Power Failure/High Temp alarms
- Lock and Key
- Battery Back Up for Controller 48hrs
- Portholes
- LED lighting

Optional Extras

- UPS Unit with Back up battery
- Preventative Maintenance program with extended warranty
- AQ Box Cloud Monitoring – with SMS and email alerts of power and temperature breaches

DIMENSIONS (HxWxD)	2060 x 595 x 633
GROSS WEIGHT	140 Kg
GROSS CAPACITY	256Lt.
TEMPERATURE RANGE	-60° to -86°C
NUMBER OF DRAWERS/SHELVES	5
POWER CONSUMPTION PER 24Hr	10.3 kWh
REFRIGERANT	Nature R

Local distributor:

Team Medical Supplies
Phone 1300 22 44 50
info@teammed.com.au
www.teammed.com.au

Medisafe
biomedical

